

Synspar Farbzuordnungstabelle

Farbe	Opaker	Dentin	Schneide	Opal Schneide
A1	O-A1	A1	Neutral	White
A2	O-A2	A2	Light	Light
A3	O-A3	A3	Grey	Clear
A3.5	O-A3.5	A3.5	Grey	Clear
A4	O-A4	A4	Grey	Grey
B1	O-B1	B1	Neutral	White
B2	O-B2	B2	Light	Light
B3	O-B3	B3	Grey	Clear
B4	O-B4	B4	Grey	Clear
C1	O-C1	C1	Light	Grey
C2	O-C2	C2	Grey	Clear
C3	O-C3	C3	Grey	Clear
C4	O-C4	C4	Grey	Grey
D2	O-D2	D2	Grey	Clear
D3	O-D3	D3	Grey	Clear
D4	O-D4	D4	Grey	Grey